

AT-Nr. / AT-No.	Ref-Nr./ ID-No. *	Hersteller / Manufacturer	Testname / Test name	Zielantigen / target antigen	Sensitivität / Sensitivity			Gesamt- Sensitivität / total sensitivity
					Cq ≤25	Cq 25-30	Cq ≥30	
AT019/20	COVAG025	GenBody Inc.	GenBody COVID-19 Ag	N	94,4%	26,1%	0,0%	46,0%
AT1100/21	CoVSLFA-20	Genobio Pharmaceutical Co., Ltd.	Virusee® SARS-CoV-2 Antigen Rapid Test (Colloidal Gold) (Saliva/Swab)	N	100,0%	95,0%	40,0%	86,0%
AT111/20	52026075	Genrui Biotech Inc.	Genrui SARS-CoV-2 Antigen Test Kit (Colloidal Gold)	N	94,1%	56,5%	0,0%	58,0%
AT161/20	P2004	GenSure Biotech Inc.	DIA-COVID® COVID-19 Ag Rapid Test Kit	N	94,1%	13,0%	0,0%	38,0%
AT166/20	CG2061	Getein Biotech, Inc.	One Step Test for SARS-CoV-2 Antigen (Colloidal Gold)	N	100,0%	82,6%	0,0%	72,0%
AT755/21	600008	Glallergen CO., LTD.	Novel Corona Virus (2019-nCoV) Ag Rapid Test Kit	N	100,0%	100,0%	60,0%	92,0%
AT533/21	CG123005	Goldsite Diagnostics Inc.	Goldsite COVID-19 SARS-CoV-2 Antigen Kit (Colloidal Gold)	N	100,0%	35,0%	0,0%	54,0%
AT018/20	n.a.	Green Cross Medical Science Corp. (Weko Pharma GmbH)	Genedia W Covid-19 Ag	N	83,3%	8,7%	0,0%	34,0%
AT542/20	n.a.	Guangdong Hecin Scientific, Inc.	2019-nCoV Antigen Test Kit (colloidal gold method)	N	82,4%	13,0%	0,0%	34,0%
AT731/21	LS-C-T-009	Guangdong Longsee Biomedical Co., Inc.	Longsee 2019-nCoV Ag Rapid Detection Kit (Immuno-Chromatography)	N	100,0%	100,0%	100,0%	100,0%
AT086/20	BE0040	Guangdong Wesail Biotech Co., Ltd.	COVID-19 Ag Test Kit	N	100,0%	52,2%	11,1%	62,0%
AT582/21	n.a.	Guangzhou Tebsun Bio-Tech Development Co., Ltd.	Tebsun 2019-nCoV Antigen Test Kit	N	95,0%	65,0%	0,0%	64,0%
AT125/20	W196	Guangzhou Wondfo Biotech Co. Ltd	Wondfo SARS-CoV-2 Antigen Test (Lateral Flow Method)	N	88,2%	0,0%	0,0%	30,0%
AT276/21	ICOV-802	Hangzhou AllTest Biotech Co., Ltd.	Covid-19 Antigen Rapid Test (Oral fluid)	N	100,0%	55,0%	0,0%	62,0%
AT766/21	INCP-502-N	Hangzhou AllTest Biotech Co., Ltd.	AllTest SARS-CoV-2 Antigen Rapid Test (Nasal Swab)	N	90,0%	20,0%	0,0%	44,0%
AT568/21	302281	Hangzhou Biotest Biotech Co., Ltd	Sienna COVID-19 Antigen-Schnelltestkassette (Nasenabstrich)	N	100,0%	95,0%	60,0%	90,0%
AT825/21	n.a.	Hangzhou Careomedic Tech Co., Ltd.	SARS-CoV-2-Antigen-Schnelltest	N	75,0%	0,0%	0,0%	30,0%
AT079/20	2020110090	Hangzhou Clongene Biotech Co., Ltd.	Clongene COVID-19 Antigen Rapid Test	N	94,4%	34,8%	0,0%	50,0%
AT816/21	n.a.	Hangzhou DIAN Biotechnology Co., Ltd.	GENEDIAN COVID-19 Antigen Test Cassette	N	90,0%	5,0%	0,0%	38,0%
AT415/21	n.a.	Hangzhou Funworld Biotech. Co., Ltd	SARS-CoV-2-Antigen-Schnelltestgerät	N	100,0%	75,0%	0,0%	70,0%
AT818/21	P211128	Hangzhou Genesis Biodetection & Biocontrol Co., Ltd.	KaiBiLi COVID-19 Antigen Pro	N	100,0%	90,0%	0,0%	76,0%
AT228/20	YMC050G	Hangzhou Immuno Biotech Co., Ltd.	IMMUNOBIO SARS-CoV-2 Antigen-Schnelltest (COVID-19 Ag)	N	88,2%	13,0%	0,0%	36,0%

Stand 14.12.2021